

Product datasheet for **RC215375L3V**

LMX1A (NM_177398) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LMX1A (NM_177398) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LMX1A
Synonyms:	DFNA7; LMX1; LMX1.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177398
ORF Size:	1146 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215375).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177398.2 , NP_796372.1
RefSeq Size:	3382 bp
RefSeq ORF:	1149 bp
Locus ID:	4009
UniProt ID:	Q8TE12
Cytogenetics:	1q23.3
Protein Families:	ES Cell Differentiation/IPS, Transcription Factors
MW:	42.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a homeodomain and LIM-domain containing protein. The encoded protein is a transcription factor that acts as a positive regulator of insulin gene transcription. This gene also plays a role in the development of dopamine producing neurons during embryogenesis. Mutations in this gene are associated with an increased risk of developing Parkinson's disease. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]