

## Product datasheet for **RC215385L3V**

### WDR27 (NM\_182552) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WDR27 (NM_182552) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDR27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_182552
ORF Size:	2571 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215385).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182552.2</a>
RefSeq Size:	3193 bp
RefSeq ORF:	2688 bp
Locus ID:	253769
UniProt ID:	<a href="#">A2RRH5</a>
Cytogenetics:	6q27
MW:	93.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein with multiple WD repeats. Proteins with these repeats may form scaffolds for protein-protein interaction and play key roles in cell signalling. Alternative splicing results in multiple transcript variants, but the full-length structure of some of these variants cannot be determined. [provided by RefSeq, Nov 2015]