

Product datasheet for RC215396L1

EEA1 (NM_003566) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EEA1 (NM_003566) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	EEA1
Synonyms:	MST105; MSTP105; ZFYVE2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215396).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_003566
ORF Size:	4233 bp




[View online >](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003566.2
RefSeq Size:	4874 bp
RefSeq ORF:	4236 bp
Locus ID:	8411
UniProt ID:	Q15075
Cytogenetics:	12q22
Domains:	FYVE, M, zf-C2H2
Protein Pathways:	Endocytosis
MW:	162.3 kDa
Gene Summary:	Binds phospholipid vesicles containing phosphatidylinositol 3-phosphate and participates in endosomal trafficking.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC215396L1



Double digestion of RC215396L1 using SgfI and MluI