

Product datasheet for **RC215427A1V**

Human CCR8 (NM_005201) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CCR8 (NM_005201) AAV Particle
Tag:	Myc-DDK
Symbol:	CCR8
Synonyms:	CC-CKR-8; CCR-8; CDw198; CKRL1; CMKBR8; CMKBRL2; CY6; GPRCY6; TER1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC215427 representing NM_005201 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATTATACACTTGACCTCAGTGTGACAACAGTGACCGACTACTACTACCCTGATATCTTCTCAAGCC
CCTGTGATGCGGAACCTATTTCAGACAAATGGCAAGTTGCTCCTTGCTGTCTTTATTGCCTCCTGTTTGT
ATTCAGTCTTCTGGGAAACAGCTGGTCATCCTGGTCCTTGCTGCTGCAAGAAGCTGAGGAGCATCACA
GATGTATACCTCTTGAACCTGGCCCTGTCTGACCTGCTTTTGTCTTCTCCTTCCCCTTCAGACCTACT
ATCTGCTGGACCAAGTGGTGTGGGACTGTAATGTGCAAAGTGGTGTCTGGCTTTTATTACATTGGCTT
CTACAGCAGCATGTTTTTCATCACCTCATGAGTGTGGACAGGTACCTGGCTGTTGTCCATGCCGTGTAT
GCCCTAAAGGTGAGGACGATCAGGATGGGCACAACGCTGTGCCTGGCAGTATGGCTAACCGCCATTATGG
CTACCATCCCATTGCTAGTGTATACCAAGTGGCCTCTGAAGATGGTGTCTACAGTGTATTATTTTCA
CAATCAACAGACTTTGAAGTGAAGATCTTCACCAACTTCAAAATGAACATTTTAGGCTTGTGATCCCA
TTCACCATCTTTATGTTCTGCTACATTAATAATCCTGCACAGCTGAAGAGGTGTCAAACCAACAAGA
CCAAGGCCATCAGGTTGGTGTCTATTGTGCTCATTGCATCTTTACTTTTCTGGGTCCCATTCAACGTGGT
TCTTTTCTCACTTCTTGCACAGTATGCACATCTTGGATGGATGTAGCATAAGCCAACAGCTGACTTAT
GCCACCATGTACAGAAATCATTTCTTTACTCACTGCTGTGTGAACCTGTTATCTATGCTTTTGTG
GGGAGAAAGTTCAAGAAACACCTCTCAGAAATATTTTCAAGAAAGTTGCAGCCAAATCTTCAACTACCTAGG
AAGACAAATGCCTAGGGAGAGCTGTGAAAGTCATCATCTGCCAGCAGCACTCTCCCGTTCTCCAGC
GTAGACTACATTTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC215427 representing NM_005201 Red=Cloning site Green=Tags(s) MDYTLDLSTTTVDYYPDIFFSSPCDAELIQTNGKLLAVFYCLLFVFSLLGNSLVILVLVVCCKLRISIT DYVLLNLALSDLLFVFSFPFQTYLLDQWVFGTVMCKVVSFGFYIGFYSSMFFITLMSVDRLAVVHAVY ALKVRTIRMGTTLCLAVWLTAIMATIPLLVFYQVASEDGLQCYSFYNQTLKWKIFTNFKMNIILGLLIP FTIFMFCYIKILHQLKRCQNHNTKAIKRLVLIIVIASLLFWVPFNVVLFSLHSMHILDGCSISQQLTY ATHVTEIISFTHCCVNPVIYAFVGEKFKKHLSEIFQKSCSQIFNYLGRQMPRESCEKSSSCQHQSSRSSS VDYIL* TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005201
ORF Size:	1065 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005201.2
RefSeq Size:	2370 bp
RefSeq ORF:	1068 bp
Locus ID:	1237
UniProt ID:	P51685
Cytogenetics:	3p22.1
MW:	40.7 kDa