

Product datasheet for **RC215459L4V**

BATF2 (NM_138456) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BATF2 (NM_138456) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BATF2
Synonyms:	SARI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138456
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215459).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138456.3
RefSeq Size:	2142 bp
RefSeq ORF:	825 bp
Locus ID:	116071
UniProt ID:	Q8N1L9
Cytogenetics:	11q13.1
Protein Families:	Transcription Factors
MW:	29.2 kDa



[View online »](#)

Gene Summary:

AP-1 family transcription factor that controls the differentiation of lineage-specific cells in the immune system. Following infection, participates in the differentiation of CD8(+) thymic conventional dendritic cells in the immune system. Acts via the formation of a heterodimer with JUN family proteins that recognizes and binds DNA sequence 5'-TGA[CG]TCA-3' and regulates expression of target genes (By similarity). Selectively suppresses CCN1 transcription and hence blocks the downstream cell proliferation signals produced by CCN1 and inhibits CCN1-induced anchorage-independent growth and invasion in several cancer types, such as breast cancer, malignant glioma and metastatic melanoma. Possibly acts by interfering with AP-1 binding to CCN1 promoter.[UniProtKB/Swiss-Prot Function]