

Product datasheet for **RC215460A1V**

Human LARP6 (NM_197958) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LARP6 (NM_197958) AAV Particle
Tag: Myc-DDK
Symbol: LARP6
Synonyms: ACHN
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC215460 representing NM_197958
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCAGTCCGGCGGGGAGGCTCGGCCGGGCCAAGACGGCGGTGCAGATCCGCGTCGCCATCCAGG
 AGGCCGAGGACGTGGACGAGTTGGAGGACGAGGAGGAGGGGCGGAGACTCGGGGCGCCGGGGACCCGGC
 CCGGTACCTCAGCCCCGGCTGGGGCAGCGCAGCGAGGAGGAGCCGAGCCGGGGCACAGGAATAGGAGC
 TCGGTGAATTCAAGGACCATGCTTGCTTCGTTTCATTGTATCTTCAGCACCTAGCACAGCACCCAGCACA

A**CGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215460 representing NM_197958
 Red=Cloning site Green=Tags(s)

MAQSGGEARPGPKTAVQIRVAIQEAEDVDELEDEEEGAETRGAGDPARYLSPGWGSASEEPPSRGHRNRS
 SVNSRTMLASFIVSSAPSTAPST

TRTRPLE**QKLISEEDLA**ANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_197958
ORF Size: 279 bp
Buffer: PBS with 0.001% Pluronic F68


[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_197958.1</u>
RefSeq Size:	577 bp
RefSeq ORF:	282 bp
Locus ID:	55323
UniProt ID:	<u>Q9BRS8</u>
Cytogenetics:	15q23
MW:	9.6 kDa