

Product datasheet for **RC215474L4V**

DAPP1 (NM_014395) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DAPP1 (NM_014395) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DAPP1
Synonyms:	BAM32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014395
ORF Size:	840 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215474).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014395.2
RefSeq Size:	2953 bp
RefSeq ORF:	843 bp
Locus ID:	27071
UniProt ID:	Q9UN19
Cytogenetics:	4q23
Domains:	SH2, PH
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)

Protein Pathways: B cell receptor signaling pathway

MW: 32.6 kDa

Gene Summary: May act as a B-cell-associated adapter that regulates B-cell antigen receptor (BCR)-signaling downstream of PI3K.[UniProtKB/Swiss-Prot Function]