

## Product datasheet for **RC215475L4V**

### FAM46C (TENT5C) (NM\_017709) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM46C (TENT5C) (NM_017709) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TENT5C
Synonyms:	FAM46C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017709
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215475).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017709.2</a>
RefSeq Size:	5720 bp
RefSeq ORF:	1176 bp
Locus ID:	54855
UniProt ID:	<a href="#">Q5VWP2</a>
Cytogenetics:	1p12
MW:	44.9 kDa



[View online »](#)

**Gene Summary:**

Nucleotidyltransferase that act as a non-canonical poly(A) RNA polymerase which enhances mRNA stability and gene expression. Mainly targets mRNAs encoding endoplasmic reticulum-targeted protein and may be involved in induction of cell death.[UniProtKB/Swiss-Prot Function]