

## Product datasheet for **RC215491L4V**

### IL26 (NM\_018402) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL26 (NM_018402) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IL26
Synonyms:	AK155; IL-26
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018402
ORF Size:	513 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215491).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018402.1</a>
RefSeq Size:	1047 bp
RefSeq ORF:	516 bp
Locus ID:	55801
UniProt ID:	<a href="#">Q9NPH9</a>
Cytogenetics:	12q15
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)

**MW:** 19.8 kDa

**Gene Summary:** This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri. [provided by RefSeq, Jul 2008]