

Product datasheet for **RC215499**

ZFP200 (ZNF200) (NM_198088) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZFP200 (ZNF200) (NM_198088) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZFP200
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215499 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGCTGCAAAAGTGGTTCCTATGCCCCAAAGCCAAAGCAGTCCTTTATACTGAGAGTTCGCCAG
ACTCCAAGCTGGGCCAAGACCTACTTCGAGATGCCACTAACGGGCCAAGACCATCCACCAGCTAGTGCT
GGAGCACTTCCTCACCTCTTGCCCAAGCCAAGCCTGGTCCAGCCAGTCAGAAAGTCAAGGAGACCTTG
GTTATTATGAAAGATGTGAGCTCAAGCCTTCAGAACAGAGTGCATCCTCGTCCCTTGGTGAAGCTTCTGC
CCAAAGGAGTCCAAAAGGAACAAGAGACAGTGTCTCTGTATTTGAAAGCTAACCCGAGGAGCTGGTGG
CTTTGAGGATTTGAATGTATTTCACTGCCAGGAAGAATGTGTGAGCTTGGATCCTACTCAACAACCTCACG
TCAGAGAAGGAAGATGACAGCAGTGTGGGAAATGATGTTACTGGCAGTCAATGGCAGTAATCCTGAAG
GTGAAGATCCTGAGAGGGAACTGTAGAAAATGAAGATTATAGAGAAAAGTCTTCAGATGATGATGAAAT
GGATTCTTCTGGTCTCTCAGCAGCCTCCCGATAACCAGGAAAAGGAACGACTAAATACATCCATTCCA
CAAAAAGGAAAATGAGAAAATCTGTTAGTTACCATTGAGAATGATACTCCTCTAGAGGAACCTCAAAT
ATGTAGACATCAGTATTATGTCCTTACTCGAAATCGGAGGACAAGGAGATGGTACACTTGTCCACTGTG
TGGGAAACAGTTAATGAAAGTTCTTACCTCATTTCCACCAGAGACCCACACTGGAGAAAAACCCAT
GACTGTAATCACTGTGGGAAAAGCTTCAATCATAAAACAACCTCAATAAACATGAGCGAATTCATACAG
GAGAGAAAACCTTATTCTGTTCTCAGTGTGGAAAAAATTCGGTCAAGAAATTCATCGGAGTCGTCATGA
AGGAATCCATATAAGGGAGAAGATATTTAAGTGTCCAGAATGTGGGAAAACCTTCCCAAAGAATGAGGAG
TTTGTGCTTATCTGCAGAGTCATGAGGCTGAGAGACCATATGGTTGCAAAAAATGTGGGAGAAGATTTG
GTCGGCTGTCAAACGTACCCGGCATGAGAAAACCCACTCAGCCTGTAAGACCCGAAAGCAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC215499 protein sequence
Red=Cloning site Green=Tags(s)

MMAAKVVPMPKPKQSFILRVPPDSKLGQDLLRDATNGPKTIHQVLVLEHFLTFLPKPSLVQPSQKVKETL
 VIMKDVSSSLQNRVHPRPLVKLLPKGVQKEQETVSLYKANKPEELVVFEDLNVFHCQEECVSLDPTQQLT
 SEKEDDSSVGGEMMLLAVNGSNPEGEDPEREPVENEDYREKSSDDDEMDSSLVSVQPPDNQEKERLNTSIP
 QKRKMRNLLVTIENDTPEEELSKYVDISIIALTRNRRTRRWYTCPLCGKQFNESYLISHQRTHTGKPY
 DCNHCGKSFNHKTNLNKHERIHTGKPYSCSQCGKNFRQNSHRSRHEGIHIREKIFKCPECGKTFPKNEE
 FVLHLQSHEAERPYGCKKCGRRFGRLSNCTRHEKTHSACKTRKQK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6379_c03.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_198088

ORF Size: 1185 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198088.3](#)

RefSeq Size: 3306 bp

RefSeq ORF: 1188 bp

Locus ID: 7752

UniProt ID: [P98182](#)


Cytogenetics: 16p13.3

Protein Families: Transcription Factors


MW: 45.5 kDa

Gene Summary: Could have a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC215499



Western blot validation of overexpression lysate (Cat# [LY405075]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215499 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).