

Product datasheet for **RC215539**

CT47A5 (NM_001080142) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CT47A5 (NM_001080142) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CT47A5
Synonyms:	CT47.5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215539 representing NM_001080142 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTCTGCCACAGGGGACCGACACCCGACCCAAGGGGACCAGGAGGCCCGGTAAAGCCAGGAGGGAGCAC
AGGCCGAGGCGGCCGAGCTGGTAACCAGGAGGGCGCGACTCCGGCCCCGACAGCAGCGACGTGGTGCC
TGCGGCCGAGGTGGTCGGAGTCGCAGGGCCCGTGAAGGCCTCGGGAGGAGGAGGGTGAGCAGGCGGCA
GGCCTGGCCGAGTCCCCGGGGCGGGAGCGCCGAGGAGGACTCAGATATCGGGCCCGCAGCGAGGAAG
AGGAGGAGGAAGAGGGGAACGAGGCGGCCAACTTCGACTTGGCGGTGGTCGCCCGTCGCTACCCGGCGTC
GGGCAATCACTTCGTGCTCCTGGACATGGTCCACTCCCTTCTCCACCGCCTCTCTACAACGACCACATC
CTCATAGAGAACCGTCAACTCAGCCGCCTGATGGTGGGGCCACACGCTGCTGCGCGCAACCCCTGGGGCA
ACCTCCCCCGCTGCTGCTGCCCCAGAGGCTGGGTGCAGGGGCCGAGCCCGGGCGGGCGAGGGCCTGGG
CCTGATCCAGGAGGCCGATCGGTCCCAGAGCCTGCAGTGCCAGCTGACCTGGCCGAGATGGCCAGGGAG
CCCGCGGAGGAGGCCGAGAGGAGAAGCTCTCAGAGGAGGCCACAGAGGAACGAGCGCAGAGGAACCGG
CCACAGAAGAACCGACCGCACAGGAGGCCACGGCCCCAGAGGAAGTCACTAAATCTCAGCCCGAAAAGTG
GGATGAAGAGGCCCAAGATGCTGCAGGCGAGGAAGAGAAAACAAGAAAAAGAGAAGGATGCGGAAAAC
AAGGTGAAGAACTCCAAGGGACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC215539 representing NM_001080142
 Red=Cloning site Green=Tags(s)


MSATGDRHPTQGDQEAPVSQEGAQAEAAGAGNQEGGDSGPDSDDVVPAAEVVGVAGPVEGLGEEEGEQAA
 GLAAVPRGGSAAEDSDIGPATEEEEEEGNEAANFDLAVVARRYPASGIHFVLLDMVHSLHRLSHNDHI
 LIENRQLSRLMVGPHAAARNPWGNLPPLLLPQRLGAGAAARAGEGLGLIQEAASVPEPAVPADLAEMARE
 PAEEAAEEKLSEEATEEPDAEEPATEEPTAQEATAPEEVTKSQPEKWDEEAQDAAGEEEEKQEKEKDAEN
 KVKNSKGT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1469_d02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001080142

ORF Size: 864 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001080142.1, NP_001073611.1

RefSeq Size: 1294 bp

RefSeq ORF: 867 bp


Locus ID: 728072

UniProt ID: Q5JQC4


Cytogenetics: Xq24

MW: 29.9 kDa

Product images:



Circular map for RC215539



Western blot validation of overexpression lysate (Cat# [LY421598]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215539 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).