

## Product datasheet for **RC215561L4V**

### **C1orf86 (FAAP20) (NM\_182533) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C1orf86 (FAAP20) (NM_182533) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1orf86
Synonyms:	C1orf86; FP7162
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182533
ORF Size:	540 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182533.2</a>
RefSeq Size:	745 bp
RefSeq ORF:	543 bp
Locus ID:	199990
UniProt ID:	<a href="#">Q6NZ36</a>
Cytogenetics:	1p36.33
MW:	19.6 kDa



[View online »](#)

**Gene Summary:**

Component of the Fanconi anemia (FA) complex required to recruit the FA complex to DNA interstrand cross-links (ICLs) and promote ICLs repair. Following DNA damage recognizes and binds 'Lys-63'-linked ubiquitin generated by RNF8 at ICLs and recruits other components of the FA complex. Promotes translesion synthesis via interaction with REV1.[UniProtKB/Swiss-Prot Function]