

Product datasheet for **RC215584L3V**

GFOD1 (NM_018988) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GFOD1 (NM_018988) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GFOD1
Synonyms:	ADG-90; C6orf114
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018988
ORF Size:	1170 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215584).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018988.1
RefSeq Size:	3033 bp
RefSeq ORF:	1173 bp
Locus ID:	54438
UniProt ID:	Q9NXC2
Cytogenetics:	6p23
Domains:	GFO_IDH_MocA, GFO_IDH_MocA_C
MW:	43 kDa


[View online »](#)