

## Product datasheet for **RC215592L2V**

### **TMS1 (PYCARD) (NM\_013258) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TMS1 (PYCARD) (NM_013258) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMS1
Synonyms:	ASC; CARD5; TMS; TMS-1; TMS1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_013258
ORF Size:	585 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215592).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013258.3</a>
RefSeq Size:	936 bp
RefSeq ORF:	588 bp
Locus ID:	29108
UniProt ID:	<a href="#">Q9ULZ3</a>
Cytogenetics:	16p11.2
Domains:	PAAD_DAPIN
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Cytosolic DNA-sensing pathway, NOD-like receptor signaling pathway

**MW:** 21.4 kDa

**Gene Summary:** This gene encodes an adaptor protein that is composed of two protein-protein interaction domains: a N-terminal PYRIN-PAAD-DAPIN domain (PYD) and a C-terminal caspase-recruitment domain (CARD). The PYD and CARD domains are members of the six-helix bundle death domain-fold superfamily that mediates assembly of large signaling complexes in the inflammatory and apoptotic signaling pathways via the activation of caspase. In normal cells, this protein is localized to the cytoplasm; however, in cells undergoing apoptosis, it forms ball-like aggregates near the nuclear periphery. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]