

Product datasheet for **RC215624L4V**

LUZP4 (NM_016383) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LUZP4 (NM_016383) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LUZP4
Synonyms:	CT-8; CT-28; CT28; HOM-TES-85
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016383
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215624).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016383.2
RefSeq Size:	1754 bp
RefSeq ORF:	942 bp
Locus ID:	51213
UniProt ID:	Q9P127
Cytogenetics:	Xq23
Protein Families:	Druggable Genome
MW:	35.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a leucine-zipper protein that was first defined as a cancer testis antigens. The encoded protein is an RNA binding protein that interacts with the mRNA export receptor nuclear RNA export factor 2. Alternate splicing results in multiple transcript variants.
[provided by RefSeq, Jan 2016]