

Product datasheet for RC215628A1V

Human EDA (NM_001005610) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human EDA (NM_001005610) AAV Particle
Tag: Myc-DDK
Symbol: EDA
Synonyms: ECTD1; ED1; ED1-A1; ED1-A2; EDA-A1; EDA-A2; EDA1; EDA2; HED; HED1; ODT1; STHAGX1; TNLG7C; XHED; XLHED
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC215628 representing NM_001005610
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGCTACCCGGAGGTGGAGCGCAGGGAACCTCTGCCTGCAGCAGCGCCGCGGGAGCGAGGGAGCCAGG
GCTGCGGGTGTGGCGGGGCCCTGCCCGGGCGGGCGAAGGGAACAGCTGCCTGCTCTTCCTGGGTTTCTT
TGGCCTCTCGCTGGCCCTCCACCTGCTGACGTTGTGCTGCTACCTAGAGTTGCGCTCGGAGTTGCGGCGG
GAACGTGGAGCCGAGTCCCGCCTTGGCGGCTCGGGCACCCTGGCACCTCTGGCACCTAAGCAGCCTCG
GTGGCCTCGACCCTGACAGCCCCATCACCAGTCACCTTGGGCAGCCGTACCTAAGCAGCAGCCATTGGA
ACCGGGAGAAGCCGACTCCACTCTGACTCCAGGACGGGCACCAGATGGCCCTATTGAATTTCTTCTTC
CCTGATGAAAAGCCATACTCTGAAGAAGAAAGTAGGCGTGTTCGCCGAATAAAAGAAGCAAAAGCAATG
AAGGAGCAGATGGTAAGTCTACTCAGTTGATCCTTTATCACTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215628 representing NM_001005610
Red=Cloning site Green=Tags(s)

MGYPEVERRELLPAAAPRERGSQGC GGAPARAGEGNSCLLFLGFFGLSLALHLLTCCYLELRSELRR
ERGAESRLGSGTPGTSGLSSLGLDPDSPITSHLQGPSPKQPLEPGEAALHSDSQDGHQMALLNFFP
PDEKPYSEESRRVRRNKRKSNEGADGKSTQLILYHF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001005610
ORF Size:	534 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005610.1</u> , <u>NP_001005610.1</u>
RefSeq Size:	1089 bp
RefSeq ORF:	408 bp
Locus ID:	1896
UniProt ID:	<u>Q92838</u>
Cytogenetics:	Xq13.1
MW:	19 kDa