

Product datasheet for **RC215630L4V**

DNAJC9 (NM_015190) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC9 (NM_015190) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC9
Synonyms:	HDJC9; JDD1; SB73
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015190
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215630).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015190.3 , NP_056005.1
RefSeq Size:	1450 bp
RefSeq ORF:	783 bp
Locus ID:	23234
UniProt ID:	Q8WXX5
Cytogenetics:	10q22.2
Domains:	Dnaj
MW:	30.4 kDa



[View online »](#)

Gene Summary:

May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8.
[UniProtKB/Swiss-Prot Function]