

## Product datasheet for **RC215652L3V**

### ASTN2 (NM\_198188) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ASTN2 (NM_198188) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ASTN2
Synonyms:	bA67K19.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198188
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215652).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198188.1</a> , <a href="#">NP_937831.1</a>
RefSeq Size:	1772 bp
RefSeq ORF:	1188 bp
Locus ID:	23245
UniProt ID:	<a href="#">O75129</a>
Cytogenetics:	9q33.1
Protein Families:	Transmembrane
MW:	44.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that is expressed in the brain and may function in neuronal migration, based on functional studies of the related astrotactin 1 gene in human and mouse. A deletion at this locus has been associated with schizophrenia. Multiple transcript variants encoding different proteins have been found for this locus. [provided by RefSeq, May 2010]