

Product datasheet for **RC215665L4V**

GPR161 (NM_153832) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | GPR161 (NM_153832) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | GPR161 |
| Synonyms: | RE2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_153832 |
| ORF Size: | 1587 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC215665). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_153832.2 |
| RefSeq Size: | 8222 bp |
| RefSeq ORF: | 1590 bp |
| Locus ID: | 23432 |
| UniProt ID: | Q8N6U8 |
| Cytogenetics: | 1q24.2 |
| Domains: | 7tm_1 |
| Protein Families: | Druggable Genome, GPCR, Transmembrane |



[View online »](#)

MW: 59 kDa

Gene Summary: The protein encoded by this gene is an orphan G protein-coupled receptor whose ligand is unknown. This gene is overexpressed in triple-negative breast cancer, and disruption of this gene slows the proliferation of basal breast cancer cells. Therefore, this gene is a potential drug target for triple-negative breast cancer. [provided by RefSeq, Mar 2017]