

Product datasheet for **RC215671A1V**

Human SFRS2B (SRSF8) (NM_032102) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SFRS2B (SRSF8) (NM_032102) AAV Particle
Tag:	Myc-DDK
Symbol:	SFRS2B
Synonyms:	DSM-1; SFRS2B; SRP46
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC215671 representing NM_032102 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCTGCGGCCGCCCCCTCCCGACGTGGACGGCATGATCACCCCTCAAGGTGGACAACCTGACCTACC
GCACCTCTCCCGACAGCTTGAGGCGCGTGTTCGAGAAGTACGGGCGCGTGGGCGACGTGTACATCCCGCG
GGAGCCCCACCAAGGCGCCCCGGGGCTTCGCTTTCGTCGCTTTCACGACCGGCGCGACGCCCAAGAC
GCCGAGGCCGCCATGGACGGGGCGGAGCTGGACGGACGCGAGCTGCGGGTGCAGGTGGCGCGCTATGGCC
GCCGGGACCTGCCCGCAGCCGCCAGGGAGAGCCACGCGGCAGGTCCAGAGGCGCGGCTACGGACGGCG
GAGCCGCAGCTACGGGCGCGGAGCCGCGAGCCCCAGGCGCGACACCGCAGCCGATCCCGGGGTCCAGC
TGCTCCAGGTCCCGCAGCCGATCTCGCTATAGGGTTCTCGCTATAGCCGGTCTCCCTACAGCCGATCTC
CTTACAGCCGGTCCGCTACAGCCGCTCTCCCTACAGCAGATCTCGCTACAGGAATCTCGCTACGGCGG
ATCTCACTACAGCTCATCTGGTTACAGTAAGTCTCGCTACAGCCGATATCACAGCAGCCGGTCTCACTCG
AAGTCTGGGTCTCCACTAGCTCTCGCTCTGCATCAACCTCAAATCGAGCTCTGCGCGACGATCCAAGT
CCTCTCGGTCTCCAGGTCTCGCTCGCGGTCCAGGTCTTCATCTATGACCAGGAGTCTCCCGGGTATC
CAAGAGGAAATCCAAGTCAAGGTGCGGATCCAAGAGGCCCCCCAAGTCTCTGAAGAGGAAGACAGATG
TCCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC215671 representing NM_032102 Red=Cloning site Green=Tags(s) MSCGRPPPDVDGMITLKVDNLTYRTSPDSLRRVFEKYGRVGDVYIPREPHTKAPRGFAFVRFHDRRDAQD AEAAMDGAELDGRELRVQVARYGRRDLPRSRQGEPRGRSRGGGYGRRSRSYGRRSRSPRRRHRSRSRGPS CSRSRSRSRYRGSRYSRSPYSRSPYSRSPYSRSPYSRSPYSRSPYSGSHYSSSGYSNSRYSRYSRSHS KSGSSTSSRSASTSKSSARRSKSSSVSRSRSRSSMTRSPPRVSKRKSKSRSRSKRPPKSPEEEGQM SS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032102
ORF Size:	846 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032102.3, NP_115285.1</u>
RefSeq Size:	4349 bp
RefSeq ORF:	849 bp
Locus ID:	10929
UniProt ID:	<u>Q9BRL6</u>
Cytogenetics:	11q21
MW:	32.7 kDa