

Product datasheet for **RC215679L4V**

DZIP1 (NM_014934) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DZIP1 (NM_014934) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DZIP1
Synonyms:	DZIP; DZIpt1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014934
ORF Size:	2601 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215679).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014934.2
RefSeq Size:	7438 bp
RefSeq ORF:	2547 bp
Locus ID:	22873
UniProt ID:	Q86YF9
Cytogenetics:	13q32.1
Domains:	zf-C2H2
MW:	98.66 kDa



[View online »](#)

Gene Summary:

May participate in spermatogenesis via its interaction with DAZ (PubMed:15081113). Has a role in primary cilium formation (PubMed:19852954).[UniProtKB/Swiss-Prot Function]