

Product datasheet for **RC215699L3V**

NARG1L (NAA16) (NM_018527) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NARG1L (NAA16) (NM_018527) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NARG1L
Synonyms:	NARG1L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018527
ORF Size:	933 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018527.2
RefSeq Size:	1692 bp
RefSeq ORF:	936 bp
Locus ID:	79612
UniProt ID:	Q6N069
Cytogenetics:	13q14.11
Protein Families:	Druggable Genome, Transcription Factors
MW:	36.7 kDa



[View online »](#)

Gene Summary: Auxillary subunit of the N-terminal acetyltransferase A (NatA) complex which displays alpha (N-terminal) acetyltransferase activity.[UniProtKB/Swiss-Prot Function]