

Product datasheet for **RC215709**

Androgen Receptor (AR) (NM_001011645) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Androgen Receptor (AR) (NM_001011645) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Androgen Receptor |
| Synonyms: | AIS; AR8; DHTR; HUMARA; HYSPI; KD; NR3C4; SBMA; SMAX1; TFM |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >RC215709 representing NM_001011645 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATACTCTGGCTTACAGTTTGGAGACTGCCAGGGACCATGTTTTGCCATTGACTATTACTTTCCAC
CCCAGAAGACCTGCCTGATCTGTGGAGATGAAGCTTCTGGGTGCTACTATGGAGCTCTCACATGTGGAAG
CTGCAAGTCTTCTTCAAAGAGCCGCTGAAGGGAAACAGAAGTACCTGTGCGCCAGCAGAAATGATTGC
ACTATTGATAAATCCGAAGGAAAAATGTCCATCTTGTCTTTCGGAAATGTTATGAAGCAGGGATGA
CTCTGGGAGCCCGAAGCTGAAGAACTTGGTAATCTGAAACTACAGGAGGAAGGAGAGGCTTCCAGCAC
CACCAGCCCACTGAGGAGACAACCCAGAAGCTGACAGTGTACACATTGAAGGCTATGAATGTCAGCCC
ATCTTTCTGAATGTCCTGGAAGCCATTGAGCCAGGTGTAGTGTGTGCTGGACACGACAACAACAGCCCG
ACTCCTTTGCAGCCTTGCTCTAGCCTCAATGAACTGGGAGAGAGACAGCTTGTACACGTGGTCAAGTG
GGCCAAGGCCCTTGCTGGCTTCCGCAACTTACACGTGGAGCAGCAGATGGCTGTCATTAGTACTCCTGG
ATGGGGCTCATGGTGTGGCCATGGGCTGGCGATCCTTCCCAATGTCAACTCCAGGATGCTCTACTTCCG
CCCCTGATCTGGTTTTCAATGAGTACCGCATGCACAAGTCCCGGATGTACAGCCAGTGTGCCGAATGAG
GCACCTCTCTCAAGAGTTTGGATGGCTCCAATCAGGAAATCCCTGTGCATGAAAGCACTGCTA
CTCTTCAGCATTATCCAGTGGATGGGCTGAAAAATCAAAAATTTTGGATGAACCTCGAATGAACACTACA
TCAAGGAACTCGATCGTATCATTGCATGCAAAAGAAAAATCCACATCCTGCTCAAGACGCTTCTACCA
GCTCACCAGCTCCTGGACTCCGTGCAGCCTATTGCGAGAGAGCTGCATCAGTTCACCTTTTGACCTGCTA
ATCAAGTCACACATGGTGAAGCTGGACTTTCCGGAAATGATGGCAGAGATCATCTCTGTGCAAGTGCCCA
AGATCCTTTCTGGGAAAGTCAAGCCATCTATTCCACACCCAG

ACGGTACGGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC215709 representing NM_001011645
Red=Cloning site Green=Tags(s)


MILWLHSLTARDHVLPIIDYYFPPQKTCLICGDEASGCHYGALTCGSCKVFFKRAAEGKQKYL CASRNDCTIDKFRRNKNCPSRLRKCYEAGMTLGARKLKKLGNLKLQEEGEASSTTSPTTEETQKLTVSHIEGYEQPIFLNVLEAIEPGVVCAGHDNNQPDFAALLSSLNELGERQLVHVVKWAKALPGFRNLHVDDQMAVIQYSWMGLMVFAMGWRSFTNVNSRMLYFAPDLVFNEYRMHKSRMYSQCVRMRHLSQEFGLWQITPQEF LCMKALLFSIIPVDGLKNQKFFDELRMNYIKELDRIIACKRKNPTSCSRRFYQLTKLLDSVQPIARELHQFTFDLLIKSHMVSVDPFEMMAEII SVQVPKILSGKVPIYFHTQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg2704_g08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001011645


ORF Size: 1164 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|--|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001011645.3 |
| RefSeq Size: | 1765 bp |
| RefSeq ORF: | 1167 bp |
| Locus ID: | 367 |
| UniProt ID: | P10275 |
| Cytogenetics: | Xq12 |
| Protein Families: | Druggable Genome, Nuclear Hormone Receptor, Transcription Factors |
| Protein Pathways: | Oocyte meiosis, Pathways in cancer, Prostate cancer |
| MW: | 44.5 kDa |
| Gene Summary: | The androgen receptor gene is more than 90 kb long and codes for a protein that has 3 major functional domains: the N-terminal domain, DNA-binding domain, and androgen-binding domain. The protein functions as a steroid-hormone activated transcription factor. Upon binding the hormone ligand, the receptor dissociates from accessory proteins, translocates into the nucleus, dimerizes, and then stimulates transcription of androgen responsive genes. This gene contains 2 polymorphic trinucleotide repeat segments that encode polyglutamine and polyglycine tracts in the N-terminal transactivation domain of its protein. Expansion of the polyglutamine tract from the normal 9-34 repeats to the pathogenic 38-62 repeats causes spinal bulbar muscular atrophy (SBMA, also known as Kennedy's disease). Mutations in this gene are also associated with complete androgen insensitivity (CAIS). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2017] |


Product images:



Circular map for RC215709



DNA-binding activity of Androgen Receptor (AR) was measured in OriGene over-expression lysate [LY423281] and a control lysate. Three microliters of each lysate was tested with a transcription factor binding assay utilizing AR-specific DNA sequences. The high level of activity observed in the over-expression lysate compared to the control lysate demonstrates that the expressed AR is biologically active in the lysate. Overexpression cell lysates are prepared from HEK293T cells transfected with RC215709 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Western blot validation of overexpression lysate (Cat# [LY423281]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215709 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).