

Product datasheet for RC215745

PANK2 (NM_153640) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | PANK2 (NM_153640) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | PANK2 |
| Synonyms: | C20orf48; HARP; HSS; NBIA1; PKAN |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >RC215745 representing NM_153640 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCCTGCTTTTATTCAAATGGGCAGAGATAAAAACTTCTCGAGTCTCCACACTGTCTTTTGTGCCACTG
GAGGTGGAGCGTACAAATTTGAGCAGGATTTTCTACAATAGGTGATCTCAGCTTTGCAAAGTGGATGA
ACTAGATTGCTTGATCAAAGGAATTTTATACATTGACTCAGTCGGATTCAATGGACGGTCACAGTGCTAT
TACTTTGAAAACCCTGCTGATTCTGAAAAGTGTGAGAAGTTACCATTTGATTTGAAAAATCCGTATCCTC
TGCTTCTGGTGAACATTGGCTCAGGGGTTAGCATCTTAGCAGTATATTCCAAAGATAATTACAAACGGGT
CACAGGTACTAGTCTGGAGGAGAACTTTTTTGGTCTCTGCTGCTTCTTACTGGCTGTACCCTTTT
GAAGAAGCTCTTGAAATGGCATCTCGTGGAGATAGCACCAAAGTGGATAAACTAGTACGAGATATTTATG
GAGGGGACTATGAGAGGTTTGGACTGCCAGGCTGGGCTGTGGCTTCAAGCTTTGAAACATGATGAGCAA
GGAGAAGCGAGAGGCTGTGAGTAAAGAGGACCTGGCCAGAGCGACTTTGATCACCATCACCAACAACATT
GGCTCAATAGCAAGAATGTGTGCCCTTAATGAAAACATTAACCAGGTGGTATTTGTTGAAAATTTCTTGA
GAATTAATACGATCGCCATGCGGCTTTTGGCATATGCTTTGGATTATTGGTCCAAGGGGCAGTTGAAAGC
ACTTTTTTCGGAACACGAGGGTTATTTTGGAGCTGTTGGAGCACTCCTTGAGCTGTTGAAGATCCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC215745 representing NM_153640
 Red=Cloning site Green=Tags(s)

MPAFIQMGRDKNFSSLHTVFCATGGGAYKFEQDFLTIGDLQLCKLDELDCIKGILYIDSVGFNGRSQCY
 YFENPADSEKQKLPFDLKNPYLLL VNIIGSGVSILAVYSKDNYKRVGTSLGGGTFGLCCLLTGCTTF
 EEALEMASRGDSTKVDKLV RDIYGGDYERFGLPGWAVASSFGNMSKEKREAVSKEDLARATLITITNNI
 GSIARMCALNENINQVVFVGNFLRINTIAMRLLAYALDYWSKGQLKALFSEHEGYFGAVGALLELLKIP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6706_f11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_153640

ORF Size: 837 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

| | |
|-------------------------------|---|
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_153640.3 , NP_705904.1 |
| RefSeq Size: | 1768 bp |
| RefSeq ORF: | 840 bp |
| Locus ID: | 80025 |
| UniProt ID: | Q9BZ23 |
| Cytogenetics: | 20p13 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Metabolic pathways, Pantothenate and CoA biosynthesis |
| MW: | 30.6 kDa |
| Gene Summary: | <p>This gene encodes a protein belonging to the pantothenate kinase family and is the only member of that family to be expressed in mitochondria. Pantothenate kinase is a key regulatory enzyme in the biosynthesis of coenzyme A (CoA) in bacteria and mammalian cells. It catalyzes the first committed step in the universal biosynthetic pathway leading to CoA and is itself subject to regulation through feedback inhibition by acyl CoA species. Mutations in this gene are associated with HARP syndrome and pantothenate kinase-associated neurodegeneration (PKAN), formerly Hallervorden-Spatz syndrome. Alternative splicing, involving the use of alternate first exons, results in multiple transcripts encoding different isoforms. [provided by RefSeq, Jul 2008]</p> |

Product images:



Circular map for RC215745



Western blot validation of overexpression lysate (Cat# [LY410958]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC217228] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).