

Product datasheet for RC215748

PTPN20B (PTPN20) (NM_001042357) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPN20B (PTPN20) (NM_001042357) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PTPN20
Synonyms:	bA42B19.1; bA142I17.1; CT126; PTPN20A; PTPN20B
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC215748 representing NM_001042357 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTTCACCTAGGGACTTTAGAGCAGAGCCTGTAACGATTATGAGGGAAATGACTCTGAAGCAGAAG
ACTTGAATTTTCAGGGAGACTTTGCCTTCATCAAGTCAGGAAAACACACCTAGATCAAAGGTTTTTGAAAA
TAAAGTTAATTCAGAGAAGGTAAAACCTTCTCTCGGAATTTCCACATAATGATTATGAGGATGTTTTT
GAAGAGCCTTCAGAAAAGTGGCAGTGATCCCAGCATGTGGACAGCCAGAGGCCCTTCAGAAGAGACAGGT
GGAGCAGTGAGGATGAGGAGGCTGCAGGGCCATCACAGGCTCTCTCCCCTCTACTTTCTGATACGCGCAA
AATTGTTTCTGAAGGAGAACTAGATCAGTTGGCTCAGATTCGGCCATTAATATTCAATTTTCATGAGCAG
ACAGCCATCAAGGATTGTTTAAAAATCCTTGAGGAAAAACAGCAGCGTATGATATCATGCAGGAATTTA
TGGCTTTAGAAGTTAAGAATCTGCCTGGTGAGTTCAACTCTGGGAATCAACCAAGCAACAGAGAAAAAA
TAGATACCGAGATATTCTCCATATGATTCAACACGCGTTCCTCTTGAAAAAGCAAGGACTACATCAAT
GCTAGTTATATTAGAATAGTCAATTGTGGAGAAGAGATTTTTATATCGCTACTCAAGGACCACCTGCTGA
GCACCATAGATGACTTTTGGCAAATGGTGGAAAAATAATTCAAATGTTATTGCCATGATAACCCAGAGA
GATAGAAGGTGGAATTATCAAATGCTACCATTACTGGCCATTTCTCTGAAGAAGCCATTGGAATTGAAA
CACTTCCGTGATTCTCGGAGAACTACAGATACTTCAATATTTTCATCATTCCAATGTTTCAAGTTGTGG
AGAAGTCCACGGAACTAGTCACTCTGTAAAACAGTTGCAAGTTCACCAAGTGGCCAGACCATGGCACTCC
TGCCTCAGCAGATAGCTTCATAAAATATATTGTTATGCAAGGAAGAGCCACCTTACAGGACCCATGGTT
GTTCACTGCAGTGCCGGCATAGGCCGGACAGGGGTGTTCTATGTGTGGATGTCGTGTTCTGTGCCATCG
TAAAGAAGTGTTCATTCAACATCATGGATATAGTGGCCCAAATGAGAGAAACACGTTCTGGCATGGTTCA
AACGAAGGAGCAGTATCACTTTTGTACGATATTGTGCTTGAAGTTCTTCGAAACTTCTGACTTTGGAT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC215748 representing NM_001042357
 Red=Cloning site Green=Tags(s)

MSSPRDFRAEPVNDYEGNDSEAEDLNFRETLPSSSQENTPRSKVFENKVNSEKVKLSLRNFPHNDYEDVF
 EEPSESGSDPSMWTARGPFRDRWSSSEDEEAAGPSQAL SPLSDTRKIVSEGELDQLAQIRPLIFNFHEQ
 TAIKDKLKILEEKTAAYDIMQEFMALELKNLPGEFNSGNQPSNREKNRYRDILPYDSTRVPLGKSKDYIN
 ASYIRIVNCGEEYFYIATQGPLLSTIDDFWQMVLENNSNVIAMI TREIEGGI IKCYHYWPI SLKKPLELK
 HFRVFL ENYQILQYFIIRMFQVVEKSTGTSHSVKQLQFTKWPDHGTPASADSF IKYIRYARKSHL TGPMV
 VHCSAGIGRTGVFLCVDVVFCAIVKNC SFNIMDIVAQMREQRSGMVQTKEQYHFCYDIVLEVLRLKLLTLD


SGP TRTRRL EQKLI SEEDLA ANDILDYK DDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001042357

ORF Size:	1260 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001042357.4
RefSeq Size:	2965 bp
RefSeq ORF:	1263 bp
Locus ID:	26095
UniProt ID:	Q4JDL3
Cytogenetics:	10q11.22
Protein Families:	Druggable Genome, Phosphatase
MW:	48.2 kDa
Gene Summary:	The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]