

Product datasheet for **RC215771L3V**

Septin 3 (SEPT3) (NM_145733) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Septin 3 (SEPT3) (NM_145733) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Septin 3
Synonyms:	bK250D10.3; SEP3; SEPT3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145733
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215771).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145733.2
RefSeq Size:	2469 bp
RefSeq ORF:	1077 bp
Locus ID:	55964
UniProt ID:	Q9UH03
Cytogenetics:	22q13.2
Domains:	GTP_CDC
MW:	40.5 kDa



[View online »](#)

Gene Summary:

This gene belongs to the septin family of GTPases. Members of this family are required for cytokinesis. Expression is upregulated by retinoic acid in a human teratocarcinoma cell line. The specific function of this gene has not been determined. Alternative splicing of this gene results in several transcript variants encoding different isoforms. [provided by RefSeq, May 2018]