

Product datasheet for **RC215773L4V**

PGRPL (PGLYRP2) (NM_052890) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PGRPL (PGLYRP2) (NM_052890) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PGRPL
Synonyms:	HMFT0141; PGLYRP2; PGRP-L; PGRPL; tagL; tagL-alpha; tagL-beta; TAGL-like
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052890
ORF Size:	1728 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215773).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052890.3
RefSeq Size:	1901 bp
RefSeq ORF:	1731 bp
Locus ID:	114770
UniProt ID:	Q96PD5
Cytogenetics:	19p13.12
Domains:	Ami_2, PGRP
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

MW: 62 kDa

Gene Summary: This gene encodes a peptidoglycan recognition protein, which belongs to the N-acetylmuramoyl-L-alanine amidase 2 family. This protein hydrolyzes the link between N-acetylmuramoyl residues and L-amino acid residues in bacterial cell wall glycopeptides, and thus may play a scavenger role by digesting biologically active peptidoglycan into biologically inactive fragments. [provided by RefSeq, Sep 2011]