

Product datasheet for RC215799A1V

Human DDX26 (INTS6) (NM_001039938) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DDX26 (INTS6) (NM_001039938) AAV Particle
Tag: Myc-DDK
Symbol: DDX26
Synonyms: DBI-1; DDX26; DDX26A; DICE1; HDB; INT6; Notch12
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC215799 representing NM_001039938
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCCATCTTACTGTTCTGATAGACACGTCTGCCTCTATGAACCAGCGCAGCCATCTGGGCACCACT
 ACCTGGACACGGCCAAAGGCGCGGTAGAGACCTTCATGAAGCTCCGTGCCCGGGACCTGCCAGCAGAGG
 AGACAGGTATATGCTGGTCACCTTCGAAGAGCCGCCCTATGCTATCAAGGCTGGATGAAAGAAAACCAT
 GCAACGTTTATGAATGAATTGAAAACTTCAGGCTGAAGGACTTACGACTCTTGCCAATCCCTAAGGA
 CAGCTTTTGATTATTAATTTAAATAGATTAGTAAGTGGCATAGACAACCTATGGGCAGGTAGGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215799 representing NM_001039938
 Red=Cloning site Green=Tags(s)

MPILLFLIDTSASMNQRSHLGGTYLDTAKGAVETFMKLRARDPASRGDRYMLVTFEPPYAIKAGWKENH
 ATFMNELKNLQAEGLTTLGQSLRTAFDLLNLNRLVTGIDNYGQVG

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001039938
ORF Size: 345 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001039938.1
RefSeq Size:	2494 bp
RefSeq ORF:	348 bp
Locus ID:	26512
UniProt ID:	Q9UL03
Cytogenetics:	13q14.3
MW:	12.7 kDa