

Product datasheet for **RC215816A1V**

Human FABP9 (NM_001080526) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human FABP9 (NM_001080526) AAV Particle

Tag: Myc-DDK

Symbol: FABP9

Synonyms: PERF; PERF15; T-FABP; TLBP

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC215816 representing NM_001080526
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTTGAGCCCTTCTTGGGAACCTGGAAGCTGGTCTCCAGTGAAAACCTTGAGGATTACATGAAAGAAC
TGGGAGTGAATTCGCAGCCCGGAACATGGCAGGGTAGTGAAACCGACAGTAACATTAGTGTTGATGG
GAAATGATGACCATAAGAACAGAAAGTTCTTCCAGGACACTAAGATCTCCTCAAGCTGGGGGAAGAA
TTTGATGAAACTACAGCAGACAACCGGAAAGTAAAGAGCACCATAACATTAGAGAATGGCTCAATGATTC
ACGTCCAAAAATGGCTTGGCAAAGAGACAACAATCAAAAGAAAAATTGTGGATGAAAAATGGTAGTGGA
ATGTAAATGAATAATTGTGCAGCACCAGAATCTACGAAAAGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215816 representing NM_001080526
Red=Cloning site Green=Tags(s)

MVEPFLGTWKLVSSENFEDYMKELGVNFAARNMAGLVKPTVTISVDGKMMTIRTESSFQDTKISFKLGEE
FDETTADNRKVKSTITLENGSMIHVQKWLKETTIRKIVDEKMMVVECKMNNIVSTRIYEKV

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001080526



[View online »](#)

ORF Size:	396 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001080526.1</u>
RefSeq Size:	399 bp
RefSeq ORF:	399 bp
Locus ID:	646480
UniProt ID:	<u>Q0Z7S8</u>
Cytogenetics:	8q21.13
MW:	14.9 kDa