

Product datasheet for **RC215829L3V**

CNKSR2 (NM_014927) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CNKSR2 (NM_014927) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CNKSR2
Synonyms:	CNK2; KSR2; MAGUIN; MRXSHG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014927
ORF Size:	3102 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215829).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014927.3
RefSeq Size:	5769 bp
RefSeq ORF:	3105 bp
Locus ID:	22866
UniProt ID:	Q8WXI2
Cytogenetics:	Xp22.12
Domains:	PDZ, SAM, PH
MW:	118 kDa



[View online »](#)

Gene Summary:

This gene encodes a multidomain protein that functions as a scaffold protein to mediate the mitogen-activated protein kinase pathways downstream from Ras. This gene product is induced by vitamin D and inhibits apoptosis in certain cancer cells. It may also play a role in ternary complex assembly of synaptic proteins at the postsynaptic membrane and coupling of signal transduction to membrane/cytoskeletal remodeling. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]