

Product datasheet for **RC215853L4V**

FKBP15 (NM_015258) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FKBP15 (NM_015258) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FKBP15
Synonyms:	FKBP133; KIAA0674; PPP1R76
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015258
ORF Size:	3657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215853).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015258.1
RefSeq Size:	4336 bp
RefSeq ORF:	3660 bp
Locus ID:	23307
UniProt ID:	Q5T1M5
Cytogenetics:	9q32
Protein Families:	Druggable Genome
MW:	133.4 kDa


[View online »](#)

Gene Summary:

May be involved in the cytoskeletal organization of neuronal growth cones. Seems to be inactive as a PPlase (By similarity). Involved in the transport of early endosomes at the level of transition between microfilament-based and microtubule-based movement.
[UniProtKB/Swiss-Prot Function]