

Product datasheet for RC215862

CEACAM19 (NM_020219) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CEACAM19 (NM_020219) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CEACAM19
Synonyms:	CEACM19; CEAL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215862 representing NM_020219 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGGAGATCCCATGGGGACCCAGGGCTGCTTCTCAAAGAGCCTCCTGCTCTCAGCCTCAATCCTGGTCC
TCTGGATGCTCCAAGGCTCCAGGCAGCTCTACATCCAGAAGATTCAGAGCAGCCTCAAAGAACCA
GGACCTTCTCCTGTCAGTCCAGGGTGTCCAGACACCTTCCAGGACTTCAACTGGTACCTGGGGAGGAG
ACGTACGGAGGCACGAGGCTATTACCTACATCCCTGGGATACAACGGCCTCAGAGGGATGGCAGTGCCA
TGGGACAGCGAGACATCGTGGGCTTCCCAATGGTTCCATGCTGCTGCGCCGCGCCAGCCTACAGACAG
TGGCACCTACCAAGTAGCCATTACCATCAACTCTGAATGGACTATGAAGGCCAAGACTGAGGTCCAGGTA
GCTGAAAAGAATAAGGAGCTGCCAGTACACACCTGCCACCAACGCTGGGATCCTGGCGCCACCATCA
TTGGATCTCTTGCTGCCGGGGCCCTTCTCATCAGCTGCATTGCCTATCTCCTGGTGACAAGGAACCTGGAG
GGGCCAGAGCCACAGACTGCCTGCTCCGAGGGGCCAGGGATCTCTGTCCATCTTGCTCGGCTGTATCC
CCAGTGCCTTCAGTGACGCCAGCACATGGATGGCGACCACAGAGAAGCCAGAATTGGGCCCTGCTCATG
ATGCTGGTGACAACAACATCTATGAAGTGATGCCCTCTCCAGTCTCCTGGTGTCCCCATCAGTGACAC
AAGGTCATAAACCCAGCCCGCCCTGCCACACCCACACCTGCAGGCGGAGCCAGAGAACCACCCAG
TACCAGCAGGACCTGCTAAACCCGACCCTGCCCTACTGCCAGCTGGTGCCAACTTCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RC215862 representing NM_020219
Red=Cloning site Green=Tags(s)



```
MEIPMGTQGCFSKSLLLSASILVLWMLQGSQAALYIQKIQEPPQKNQDLLL SVQGVPTDFQDFNWYL GEE
TYGGTRLFTYIPGIQRPQRDGSAMGQRDIVGFNGSMLLRRAQPTDSGT YQVAITINSEWTMKAKTEVQV
AEKNKELPSTHLPNAGILAATIIGSLAAGALLISCIAYLLVTRNWRGQSHR LPAPRQGGSLSILCSAVS
PVPSVTPSTWMATTEKPEL GPAHDAGDNNIYEVMPSPVLLVSPISDTRSINPARLP TTPPHLQAEPENHQ
YQQDLLNPDPAPYCQLVPTS
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8010_b09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_020219

ORF Size: 900 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020219.3](#), [NP_064604.2](#)

RefSeq Size: 1863 bp

RefSeq ORF: 903 bp

Locus ID: 56971


UniProt ID: [Q7Z692](#)

Cytogenetics: 19q13.31


Protein Families: Transmembrane

MW: 32.5 kDa


Product images:



Circular map for RC215862



Western blot validation of overexpression lysate (Cat# [LY412583]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215862 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified CEACAM19 protein (Cat# [TP315862]). The protein was produced from HEK293T cells transfected with CEACAM19 cDNA clone (Cat# RC215862) using MegaTran 2.0 (Cat# [TT210002]).