

Product datasheet for **RC215864L2V**

FA20A (FAM20A) (NM_017565) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FA20A (FAM20A) (NM_017565) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FA20A
Synonyms:	A11G; AIGFS; FP2747
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_017565
ORF Size:	1623 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215864).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017565.2
RefSeq Size:	2267 bp
RefSeq ORF:	1626 bp
Locus ID:	54757
UniProt ID:	Q96MK3
Cytogenetics:	17q24.2
Protein Families:	Secreted Protein, Transmembrane
MW:	61.2 kDa



[View online »](#)

Gene Summary:

This locus encodes a protein that is likely secreted and may function in hematopoiesis. A mutation at this locus has been associated with amelogenesis imperfecta and gingival hyperplasia syndrome. Alternatively spliced transcript variants have been identified. [provided by RefSeq, Aug 2011]