

Product datasheet for **RC215886**

RSPH6A (NM_030785) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RSPH6A (NM_030785) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RSPH6A
Synonyms:	RSHL1; RSP4; RSP6; RSPH4B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

```
>RC215886 representing NM_030785
Red=Cloning site Blue=ORF Green=Tags(s)
```

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGACCTGCCGCCCTACCTGAGCGCCCTGCCACGAGCCTCCGGGCCGGAGGAACTTCTAGGCCCTCCAGAGGGCGGCACAGTCGGGACCAAGCTCAGGCCCTGGCAGCGGACCCGAGGAGAGGCAGCAGATACCTCCAGACGCCACGCAAAACGCCCTGGTTGGTCACAGAGGGGACGCTGTCCCAACAGGAGAACTTGCTGATGCCCCAGGTCTTCCAGGCTGAGGAAGCCCGGCTGGTGGCATGGAGTACCCATCTGTGAACACGGGCTTCCCTCAGAGTTCCAGCCTCAGCCTTACTCTGATGAAAGCAGGATGCAGGTGCGCCGAGCTCACCACCAGCCTAATGCTGCAGCGGCTCCAGCAGGGCCAAAGCAGCCTGTTCCAGCAACTGGACCCACCTTCCAGGAGCCCCAGTCAACCCCTTGGGCCAGTTCACCTCTACCAGACAGACCAGTTCTCTGAAGGTGCCCAGCACGGCCTTACATAAGGGATGACCTTGCCCTTCAGTTCTTGCCCTCTGAGCTGGGCTTCCCACACTACAGTGCACAGGTGCTGAGCCCCAGCCTCTGGAGCTGGCCGTGCAGAACGCCAAGGCCTACCTGCTGCAGACCAGCATCAATTGCGACCTCAGCCTGTACGAGCACCTGGTGAATCTGCTGACCAAGATCCTGAACCAGCGGCCTGAGGACCCCTTGCTGTCTGGAGTCTCTGAACCGCACCCAGCAGTGGGAGTGGTTCACCCCAAGCTGGAACGCTGCGGGACGCCGAGCTCAGACCCAGTCAAGCAGCTACAAGATGGCGAGGAACAGAACGGCCTGTTCAACCGGAGTGGAGGCGGCACTGAAGTGCAGAACAGGAGATGGAGGAGGAGTGGGGAGACACCAAGCTGCCCCAATCATGGAGACTGCCTTCTACTTCGAGCAGGCGGCGCTCGGCCTGAGCTCGGACGAGAGCTTCCGCATTTTCTGGCCATGAAACAGCTGGTGGAGCAGCAGCCATCCACACCTGTCGTTCTGGGGCAAGATCCTGGGATCAAACGCAGCTACCTGGTGGCCGAGGTGGAATTCCGGGAGGGCGAGGAGGAGGCAGAGGAGGAGGAGTGGAGGAGATGACGGAAGGTGGCAGGTCATGGAGGCGCACGGCAGGAGGAGGGCGAGGAGGACGAGGAAGGCCGTGGACATCGTCCCTAAGTCCGTATGGAAGCCGCCGCCCGTGATCCCCAAGGAGGAGAGCCGCTCAGGCGCCAACAAGTACCTGTACTTTGTGTGCAACGAGCCGGGCTGCCATGGACGCGGCTGCCACGCTACTCCAGCCAGATCGTGAACGCCCCGAAAGATCAAGAAGTTCTTCACAGGCTACCTGGACACGCCAGTCGTAGCTACCCACCCTTCCCGGGCAACGAGGCCAACTACCTGCGGGCCAGATAGCCCGCATCTCGGCCACGACGAGGTAGCCCGCTGGGCTTCTACCAGTTTGTGAGGAGGAGGGCGACGAGGAGGAGGAAGGTGTGCTGGGCGCGACTCTACGAGGAGAACCCGGACTTCGAGGGCATCCCCGTGCTGGAGCTGGTCGACTCATGGCCAACTGGGTGCATCACACACAGCACATCTGCCGACAGGCGCTGCACTTGGGTGAACCCCTTGAGAAGACAGAGGAGGAGGAGGACCTGGGGGAGGAGGAAGAGAAGGCAGATGAGGGGCCAGAGGAGGTGAGCAGGAGGTTGGCCCCCACTGCTAACGCCACTTTCAGAAGATGCAGAAATCATGCACCTGGCACCCCTGACCAACCCCGCTGTCTGACAGCTCTGCCCGAGTACTAGTGCCGTTGTGCGCTCCAACCTCTGGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215886 representing NM_030785
 Red=Cloning site Green=Tags(s)


MGDLPPYPERPAQQPPGRRTSQASQRRHSRDQAQALAADPEERQQIPDAQRNAPGWSQRGSLSQQENLL
 MPQVFQAEERLGGMEYPSVNTGFPSEFQPPYSDSRMQVAELTSLMLQRLQQGSSLFQQLDPTFQE
 PPVNPGLGQFNLQYTDQFSEGAQHGPYIRDDPALQFLPSELGFPHYSAQVPEPEPELEAVQNAKAYLLQTS
 INCDSLSEYHLVNLTKILNQRPEPLSVLESLNRTTQWEWFHPKLDLTDDEMPQPTYKMAEKQKALFT
 RSGGGTEGEQEMEEVGETPVPNIMETAFYFEQAGVGLSSDESFRIFLAMKQLVEQQPIHTCRFWGKILG
 IKRSYLVAEVEFREGEAEAEVEEMTEGGEVMEAHGEEEGEEDDEEKAVDIVPKSVWKPVPVPIKEESR
 SGANKYL YFVCNEPGLPWTRLPHVTPAQIVNARKIKKFFTGYLDTVPVSYPPPGNEANYLRAQIARISA
 ATQVSPLGFYQFSEEGDEEEGGAGRDSYEENPDFEGIPVLELVDSMANVWHHTQHILPQGRCTWVNPL
 QKTEEEEDLGEEEKADEGPEEVEQEVGPPLLTPLEDAEIMHLAPWTTTLSCSLCPQYSVAVVRSNLWP
 GAYAYASGKKFENIYIGWGHKYSPEFNPALPAPIQQEYPSGPEIMEMSDPTVEEEQALAAEQEQALGAT
 EEEEEEEEEEEETDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8028_a12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_030785


ORF Size: 2151 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_030785.4</u>
RefSeq Size:	2477 bp
RefSeq ORF:	2154 bp
Locus ID:	81492
UniProt ID:	<u>Q9H0K4</u>
Cytogenetics:	19q13.32
MW:	80.7 kDa
Gene Summary:	The protein encoded by this gene is similar to a sea urchin radial spoke head protein. Radial spoke protein complexes form part of the axoneme of eukaryotic flagella and are located between the axoneme's outer ring of doublet microtubules and central pair of microtubules. In Chlamydomonas, radial spoke proteins are thought to regulate the activity of dynein and the symmetry of flagellar bending patterns. This gene maps to a region of chromosome 19 that is linked to primary ciliary dyskinesia-2 (CILD2). [provided by RefSeq, Jul 2008]


Product images:



Circular map for RC215886



Western blot validation of overexpression lysate (Cat# [LY410707]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215886 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified RSPH6A protein (Cat# [TP315886]). The protein was produced from HEK293T cells transfected with RSPH6A cDNA clone (Cat# RC215886) using MegaTran 2.0 (Cat# [TT210002]).