

Product datasheet for **RC215896**

FOXK1 (NM_001037165) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: FOXK1 (NM_001037165) Human Tagged ORF Clone

Tag: Myc-DDK

Symbol: FOXK1

Synonyms: FOXK1L

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide >RC215896 representing NM_001037165
 Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

GCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGCGCCAGATCTCAAGCTTAACTAGCTAGCGGACCGAC

ATGGCCGAAGTCGGCGAGGACAGCGGGCGCCGCGCCCTGCTCGCGCTGCGCTCGGCGCCCTGCAGCCAG
 TGCTGTGCGCCGACGCCGCGCCGCGCCCTTCCCGCGGCCGACCCCGCGGCCCGCGCAGCCCA
 GCCTCCGCCCGGGCGCGCCGCGCCGCCACCGCCGCTGCCTCCGGGCGCGATCGCGGGCGGGGCTCC
 TCCGGGGCTCCTCCGGGGTATCCGGGACTCCGCGGTGCGGGGCGCGCGCCGCGCCCTGGTGGCCGCG
 CGGCCGCTCGGTACGGCAGAGCCCGGGCGCGCTGGCGCGGCTGGAGGGCCGCGAGTTCGAGTTCT
 CATGCGCCAGCCAGCGTCACCATCGGCCGAACTCGTCGCGAGGGCTCGGTGGACTTGAGCATGGGCTG
 TCCAGCTTCATCTCGCGGCCACCTGCAGCTCAGCTTCCAGGAGCCGCACTTCTACCTGCGCTGCCTCG
 GCAAGAACGGCGTCTTCGTGGACGGGGCTTCCAGAGACGCGGCGCGCCCGCCCTGCAGCTGCCAAGCA
 GTGTACCTTCCGGTTTCCAGCAGGGCATCAAGATCCAGTTCACGTGCTCTATCACAAGAAGAGGCC
 CCAGCCTCCCGCTGCGGCCACTGTACCCCGAGATCTCCCTCTGAAGATCCACATCCCGGAGCCGGACC
 TCCGGAGCATGGTCAGCCCCGTCCCCTCCCGACGGGACCATCAGTGTCCCAACTCTGCCAGCCAG
 TCCACGCGGTGCCGGCTCCTCCAGTTACCGCTTTGTGAGAAGCTGACCTCGGACCTGCAGCTGGCAGCA
 GAGTTTGAGCAAAGGCCGCGTCCGAGCAGCAGGCAGACACGCTGGAGGAGACAGCCCAAGGATGAGT
 CAAAGCCGCGTTCTCTACGCGCAGCTGATCGTGCAGGCCATCTCCTCCGCCAGGACCGGCAGCTGAC
 CCTGAGCGGGATCTACGCCACATCACAAGCATTACCCCTACTACCGGACGGCCGACAAGGCTGGCAG
 AATTCTATCCGGCACAACCTCTCTTTGAACCGTTACTTTATCAAAGTCCCACGTTCCAGGAGGAGCCTG
 GGAAGGGTCTTTTGGCGAATAGACCCTGCCTCTGAAGCCAAGCTCGTGAACAGGCATTCCGAAACG
 GAGGCAGAGGGTGTCTCTGCTTCCGACCCCTTCCGGGCTGTGCTCAAGGAGCGCTCCAGCTTCG
 CCCACACACCCCGGGTGTATGTCCTCCGCTCCGGCGGCTGCAGACCCAGAGTGCCTGTCTCGGGAGG
 GCTCCCCATTCCACACGACCCTGAGTTTGGTCCAAGTTAGCTTCTGTCCAGAGTACCGGTATTCCCA
 AAGCGCACCCGGCTCCCCGTGAGCGCCAGCCAGTGCATGGCCGTGCCTCCCCGACCGTCCAGCCTC
 GTGGCAAGCCCGTGGCTACATGCCCGCTCCATCGTAACTCACAGCAGCCCGGGCCACGCCATCC
 ACGTCGTGACGAGGCCCCACCGTACCATGGTCAGGGTGGTACCACATCTGCCAATCGGCCAACGG
 ATACATCCTACCAGCCAGGGCGCGGGGGGGTCCCATGATGCGGGGGCGCAGCCGTGCTGGACCTG
 GGCAGCGAGGCCAGAGCCTGGAGGAGAAACCCACCATTGCGTTTGCCACAATCCCGCGGCTGGTGGAG
 TCATCCAGACGGTGGCCAGCCAGATGGCCCCGGGTCCCCGGACACACGGTACCATCCTGCAGCCCGC
 CACACCGTGACCTCGGGCAGCACCCTTCCAGTCCGGGCGGTGACCCAGAACGGAAGCATGCGGTT
 CCCACGAACAGTTTAGCCGGCAACGCTTACGCCCTCACCAGCCCTTTCAGCTCCTTGGACCCAAAGCA
 GTTTCATCCGCGCGGTGGTGGTCAACCGGGTGTGCGAGGTGGGGCCAAAGGAGCCAGCAGCAGCCGTCGC
 GGCCACGGCCACCACCCAGCCACTGCCACCACCGCCTCTGCCTCCGCCTTCCACTGGAGAGCCC
 GAGGTCAAAGGTCCCGGGTGGAGGAGCCAGTGGGGCTGTAACCACACCGGCTGGAGTGATCGCAGCTG
 CCGGCCCCAGGGGCCAGGCACCGGGAG

ACCGGTACGCGCGCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC215896 representing NM_001037165
 Red=Cloning site Green=Tags(s)

MAEVGEDSGARALLALRSAPCSPVLCAAAAAFAAAAPPAPAPQPQPPGPPPPPPPLPPGAIAGAGS
 SGGSSGVSGDSAVAGAAPALVAAAAA SVRQSPGPALARLEGREFEFLMRQPSVTIGRNSQGSVDLSMGL
 SSFISRRHLQLSFQEPHFYLRCLGKNGVFDGAFQRRGAPALQLPKQCTFRFPSTAIIQFTSLYHKEEA
 PASPLRPLYPQISPLKIHIPEDLRSMVSPVPSPTGTISVPNSCPASPRGAGSSSYRFVQNVTSDLQLAA
 EFAAKAASEQQADTSGGDSKPKDESKPPFSYAQLIVQAISSAQDRQLTLSGIYAHITKHYPYRTADKGWQ
 NSIRHNL SLNRYFIKIVPRSQEEPGKGSFWRIDPASEAKLVEQAFRRRQRGVSCFRTPFGPLSSRSAPAS
 PTHPGLMSRSGGLQTPCELSREGSPIPHDPEFGSKLASVPEYRYSQSAPGSPVSAQPVIMAVPPRSSL
 VAKPVAYMPASIVTSQQPAGHA IHVVQQAPTVTMVRVTT SANSANGYILTSQGAAGGSHDAAGAAVLDL
 GSEARGLEEKPTIAFATIPAAGGVIQTVASQMAPGVPGHVTIILQPATPVTLGQHHLVRAVTVNGKHAV
 PTNSLAGNAYALTSPLQLLATQASSAPVVVTRVCEVGPKEPAAVAATATTT PATATTASASASSTGEP
 EVKRSRVEEPSGAVTTPAGVIAAAGPQPGTGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: RsrII-MluI

Cloning Scheme:



ACCN: NM_001037165

ORF Size: 2199 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037165.1](#), [NP_001032242.1](#)

RefSeq Size: 11191 bp

RefSeq ORF: 2202 bp

Locus ID: 221937


UniProt ID: [P85037](#)

Cytogenetics: 7p22.1

MW: 75.3 kDa

Gene Summary: Transcriptional regulator involved in different processes such as glucose metabolism, aerobic glycolysis, muscle cell differentiation and autophagy (By similarity). Recognizes and binds the forkhead DNA sequence motif (5'-GTAAACA-3') and can both act as a transcription activator or repressor, depending on the context (PubMed:17670796). Together with FOXK2, acts as a key regulator of metabolic reprogramming towards aerobic glycolysis, a process in which glucose is converted to lactate in the presence of oxygen (By similarity). Acts by promoting expression of enzymes for glycolysis (such as hexokinase-2 (HK2), phosphofructokinase, pyruvate kinase (PKLR) and lactate dehydrogenase), while suppressing further oxidation of pyruvate in the mitochondria by up-regulating pyruvate dehydrogenase kinases PDK1 and PDK4 (By similarity). Probably plays a role in gluconeogenesis during overnight fasting, when lactate from white adipose tissue and muscle is the main substrate (By similarity). Involved in mTORC1-mediated metabolic reprogramming: in response to mTORC1 signaling, translocates into the nucleus and regulates the expression of genes associated with glycolysis and downstream anabolic pathways, such as HIF1A, thereby regulating glucose metabolism (By similarity). Together with FOXK2, acts as a negative regulator of autophagy in skeletal muscle: in response to starvation, enters the nucleus, binds the promoters of autophagy genes and represses their expression, preventing proteolysis of skeletal muscle proteins (By similarity). Acts as a transcriptional regulator of the myogenic progenitor cell population in skeletal muscle (By similarity). Binds to the upstream enhancer region (CCAC box) of myoglobin (MB) gene, regulating the myogenic progenitor cell population (By similarity). Promotes muscle progenitor cell proliferation by repressing the transcriptional activity of FOXO4, thereby inhibiting myogenic differentiation (By similarity). Involved in remodeling processes of adult muscles that occur in response to physiological stimuli (By similarity). Required to correct temporal orchestration of molecular and cellular events necessary for muscle repair (By similarity). Represses myogenic differentiation by inhibiting MEFC activity (By similarity). Positively regulates Wnt/beta-catenin signaling by translocating DVL into the nucleus (PubMed:25805136). Reduces virus replication, probably by binding the interferon stimulated response element (ISRE) to promote antiviral gene expression (PubMed:25852164).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC215896