

## Product datasheet for **RC215897L4V**

### **BTF (BCLAF1) (NM\_001077441) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BTF (BCLAF1) (NM_001077441) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BTF
Synonyms:	bK211L9.1; BTF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001077441
ORF Size:	2241 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215897).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001077441.1</a>
RefSeq Size:	6758 bp
RefSeq ORF:	2244 bp
Locus ID:	9774
UniProt ID:	<a href="#">Q9NYF8</a>
Cytogenetics:	6q23.3
Protein Families:	Transcription Factors
MW:	85.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]