

Product datasheet for **RC215898L4V**

alpha Adducin (ADD1) (NM_014189) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	alpha Adducin (ADD1) (NM_014189) Human Tagged ORF Clone Lentiviral Particle
Symbol:	alpha Adducin
Synonyms:	ADDA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014189
ORF Size:	2304 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215898).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014189.2
RefSeq Size:	4063 bp
RefSeq ORF:	2307 bp
Locus ID:	118
UniProt ID:	P35611
Cytogenetics:	4p16.3
Domains:	Aldolase_II
Protein Families:	Druggable Genome



[View online »](#)

MW: 84.1 kDa

Gene Summary: Adducins are a family of cytoskeletal proteins encoded by three genes (alpha, beta, and gamma). Adducin acts as a heterodimer of the related alpha, beta, or gamma subunits. The protein encoded by this gene represents the alpha subunit. Alpha- and beta-adducin include a protease-resistant N-terminal region and a protease-sensitive, hydrophilic C-terminal region. Adducin binds with high affinity to Ca(2+)/calmodulin and is a substrate for protein kinases A and C. [provided by RefSeq, Aug 2017]