

Product datasheet for RC215907

UEVLD (NM_018314) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UEVLD (NM_018314) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UEVLD
Synonyms:	ATTP; UEV3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215907 representing NM_018314 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTTCGACTGCGAGGGCCTGAGACGGCTGCTTGGCAAGTACAAGTTCAGGGACCTAACTGTGGAAG
AACTACGGAATGTAATGTATTTTCCACATTTCAAATATTCATGGACACCTATGTTTTAAAGATAG
TTCTCAGAAAGACCTGCTGAATTTACTGGCACAATTCCTGTGATGATCAGGGTAATACATATAACATA
CCAATTCGTTTCTGGATTTGGATTCTCACCCTTCGCTCCCCCTATTTGCTTCTTGAAGCCAAGTCAA
ATATGGGAATCTTAGTCGAAAACATGTGGATGCTCAAGGCAGAAATATATTTGCCCTATCTCCAAAAGT
GAGCCATCCTAAATCTGTCAATTGTTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAGGAAGTCCC
ATGTATTTCTATCATCATCTGATGAGGCACGGCAGGTAGACTTGCTAGCCTATATTGCAAAAATCACTG
AAGGTGTTTCAGATACAAATCAAGAGCTGGGCAATCATGAGAATAAAACAGTCAATAAAATTAAGTGT
GGTTGGAGGTGGAGAAGTGGTATTGCCTGCACATTAGCAATTCAGCAAAGGTATTGCAGACAGGCTT
GTCCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCCACGATGGACCTTGAATCTTCAACCTTCCTAATG
TGGAGATCAGCAAAGATTTGCTGCCTCTGCTCATTCCAAGGTGGTGATCTTACAGTCAACTCTTTGGG
TAGTTCTCAGTCGTACCTTGATGTGGTACAGAGCAATGTGGATATGTTTACAGGCCCTTGTCAGGCTCTG
GGACATTATAGTCAACACAGTGTCTGCTGCTGATCTCAACCAGTGGAAATCATGACCTATGTAACAT
GGAAACTGAGTACATTTCTGCAAATCGAGTGATCGGAATTGGATGTAATCTGGATTACAGAGATTACA
GTATATTATTACAAATGTTTGAAGGCACAGACTTCAGGCAAAGAAGTATGGGTTATTGGCGAGCAAGGA
GAAGACAAAGTGCTCACATGGAGTGGCCAAAGAAGTAGTGAGTCATACCTCTCAAGTGCAGCTGTCCA
ACAGGGATATTATGATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC215907 representing NM_018314
 Red=Cloning site Green=Tags(s)


MEFDCEGLRRLGKYKFRDLTVEELRNVNFFPHFKYSMDTYVFKDSSQKDLLNFTGTIPVMYQGNTYNI
 PIRFWILDSPFPAPPICFLKPTANMGILVGKHVDAQGRIYLPYLQNWHPKSVIVGLIKEMIAKFQEELP
 MYSLSSSDEARQVDLLAYIAKITEGVSDTNSKSWANHENKTVNKITVVGGEGLIACTLAISAKGIADRL
 VLLDLSEGTGATMDLEIFNLPNVEISKDLSASAHSKVVIPTVNSLGSSQSYLDVVQSNVDMFRALVPAL
 GHYSQHSVLLVASQPVEIMTYVTWKLSTFPANRVIGIGCNLDSQRLQYIITNVLKAQTSKKEVWVIGEQQ
 EDKVLTSWGSQEEVVSHTSQVQLSNRDI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg2903_e07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_018314

ORF Size: 1137 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_018314.6](#)

RefSeq Size: 1984 bp

RefSeq ORF: 1140 bp

Locus ID: 55293

UniProt ID: [Q8IX04](#)


Cytogenetics: 11p15.1

Domains: Idh


MW: 42.1 kDa

Gene Summary: Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC215907



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UEVLD (Cat# RC215907, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UEVLD (Cat# [TA809207])(1:2000). Positive lysates [LY413157] (100ug) and [LC413157] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY413157]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215907 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).