

Product datasheet for **RC215907L3V**

UEVLD (NM_018314) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UEVLD (NM_018314) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UEVLD
Synonyms:	ATTP; UEV3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018314
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018314.2
RefSeq Size:	1984 bp
RefSeq ORF:	1140 bp
Locus ID:	55293
UniProt ID:	Q8IX04
Cytogenetics:	11p15.1
Domains:	ldh
MW:	42.1 kDa


[View online »](#)

Gene Summary: Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]