

Product datasheet for **RC215908L3V**

MYZAP (NM_152451) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MYZAP (NM_152451) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MYZAP
Synonyms:	GCOM1; Gup; MYOZAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152451
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215908).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152451.1
RefSeq Size:	2356 bp
RefSeq ORF:	1317 bp
Locus ID:	100820829
UniProt ID:	P0CAP1
Cytogenetics:	15q21.3
MW:	50.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that is abundantly expressed in cardiac tissue. The encoded protein localizes to intercalated discs in cardiomyocytes and functions as an activator of Rho-dependent serum-response factor signaling. Alternative splicing results in multiple transcript variants. Readthrough transcription also exists between this gene and the neighboring downstream gene POLR2M (polymerase (RNA) II (DNA directed) polypeptide M) and is represented with GeneID: 145781. [provided by RefSeq, Mar 2014]