

Product datasheet for **RC215910**

MAGIX (NM_024859) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAGIX (NM_024859) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MAGIX
Synonyms:	JM10; PDZX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215910 representing NM_024859 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGAGCCGCGCACAGGGGGCGCCGGAACCTAAGGGGAGCAGAGGAAGCCGCGGCCCTCCCGCTCG
CGGGCCCTAGCGCCCGGAGCTCCTGGCGCGTTGGACGCGCGCCCCCTGGCGGCGCAGCTGCGGTGCA
CGTGGCAGCGCTGGTACGCAGGGCGGGCGCCACATTGCGCTGCGCCGGAAGGAGGCTGTTAGCGTGCTG
GACTCTGCGGACATAGAGGTACAGACAGTCGCCTGCCTCATGCCACTATTGTGGATCACCGGCCCCAGC
ACCGTTGGTTAGAGACATGTAACGCACCTCCCCAATTGATCCAGGGTAAGGCACGTAGTGCTCCGAAGCC
ATCCCAGGCCTCTGGTCAATTTCTCTGTGGAGCTGGTTCGCGGTTACGCAGGCTTTGGCCTCACCTTAGGT
GGGGGCCGGGATGTAGCTGGGGACACTCCGCTGGCCGTGCGCGGGCTGCTGAAGGATGGCCAGCACAGC
GCTGTGGTCTGTTGGAGGTGGGGACCTCGTGCTCCACATCAACGGAGAGTCAACGCAGGGCCTCACCCA
TGCCCAGGCCGTGGAGCGGATCCGAGCTGGAGGCCCCAGCTCCACCTGGTTATTCGTGCGCCTCTGGAG
ACCCACCTGGCAAGCCTCGAGGGGTGGGAGAGCCCCGAAAGGAGTTGTCCCGTCATGGCCAGATCGCA
GCCCAGATCCTGGAGGGCCGGAGGTAACGGGTCTCGCAGCAGCAGCACTCCCTAGTTCAGCACCTCC
ATCCCGGACGACGCTCAAGAAGACCCGGGGCAGCCCGGAGCCTAGTCCAGAGGCGGCCCGCATGGCCCC
ACGGTTTCTCCTCTGAGCGCCGCGCTGAGGATCCCAACGACCAGATCCCGGGTTCCCGGGGCCCTGGC
TAGTGCCAGCGAGGAACGGCTCTCGCGGGCCCTAGGGGTCCGGGGGCGCAGCAGCTCGCTCAGGAGAT
GGCAGCCGGAAGCGGAGACAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC215910 representing NM_024859
 Red=Cloning site Green=Tags(s)

MEPRTGGAANPKGSRGSRGPSPLAGPSARQLLARLDARPLAARAAVDVAALVRRAGATLRLRRKEAVSVL
 DSADIEVTDSRLPHATIVDHRPQHRWLETCNAPPQLIQGKARSAPKPSQASGHFSVELVRGYAGFGLTLG
 GGRDVAGDTPLA VRGLLDGPAQRGRLEVGD LVHNGESTQGLTHAQAVIRAGGPQLHLVIRRRPLE
 THPGKPRGVGEPRKGVVPSWPD RSPDPGGPEVTGSRSSSTSLVQHPPSRTTLKTRGSPEPSPEAAADGP
 TVSPERRAEDPNDQIPGSPGPWLVPSEERLSRALGVRGAAQLAQEMAAGRRRH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8028_a08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024859

ORF Size: 1002 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024859.2](#), [NP_079135.2](#)

RefSeq Size: 2150 bp

RefSeq ORF: 1005 bp


Locus ID: 79917

UniProt ID: [Q9H6Y5](#)


Cytogenetics: Xp11.23

MW: 35 kDa


Product images:



Circular map for RC215910



Western blot validation of overexpression lysate (Cat# [LY411025]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215910 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified MAGIX protein (Cat# [TP315910]). The protein was produced from HEK293T cells transfected with MAGIX cDNA clone (Cat# RC215910) using MegaTran 2.0 (Cat# [TT210002]).