

Product datasheet for **RC215918**

NUTM2F (NM_017561) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NUTM2F (NM_017561) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NUTM2F
Synonyms:	FAM22F; NUTMF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC215918 representing NM_017561
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGCTTCAAATGGAGCATACCCAGTGCTGGGACCCGGCGTGACCGTGAACCCTGGCACCTCCCTGTCTG
 TGTTACGGCTCTGCCCTTTGCCACACCCGCTCCCGGCCAGCACACAGGCCGCCCTCGTGACTGCAGT
 GGTTCCTCCAGCCGGCCCTCTGGTGCTCTCTGCCCTCCCCAGCACCCCTCTAGTGGCAGGACAGGATGGC
 CGCGGCCCGAGTGGGGCCGGGGCTTCCAACGTCTTTGTCCAGATGAGGACAGAAGTGGGGCTGTGAAGC
 CCCCTCAGGCACAGACCTTGATCCTAACTCAGGCCCCCTCGTCTGGCAGGCTCCAGGCACCTCTGTGG
 AGGTGTACATGTGTCCACCTCCCCTACTCCTGGCAGCTGCTCTGGGGTGGCCGTTACCTCTGCCAGGTG
 GTTGGGGGCACCCAGGCTGTGAGGGAGGCTGGTCCCATGGCCTTCTCTTCCACCACCACCACCGGCTG
 CCCAGGTGGCCCCATCGTGTCCCAGGGAACGCTAGGCCATGGCCACAAGGGGCTCATGGAGAGGGCAG
 CCTGGCTCCCTCCAGGCCAAGGCCCGGCCGACGACTCCTGTAACCCAAGAGTGTCTATGAGAACTTC
 CGACTCTGGCAGCACTACAAGCCCTGGCCCGGAGGCACCTTCCCAGAGTCTGACACTGAAGCGCTTT
 CCTGCTTCTCATCCAGTTCTCCGATCCCTGGCCCGGCGGAAGCCACCATGACGCTGGAGGAGGACT
 GTGGCAGGCCATGCGGGAATGGCAGCACACGAGCAACTTTGACCGGATGATCTTCTACGAGATGGCGGAA
 AAGTTCCTGGAATTTGAGGCTGAGGAGGAGATGCAGATTCAGAAATCGCAGTGGATGAAGGGGCCCGAGA
 GCCTGCCTCCTCCAGCCCCGCCGAGGCTTGAACCTCGAGGACCCCTGCCCTGAGGTGGTCAAGCAGCC
 AGTGTACCTTCCAGCAAGGATGGCCCCAAGGCCCGACTGCCTGCCTGCCACCACCAGGCCCGAGAGG
 CCAGCGGAGACCAAGGCCACCTGCCACCACCAGGCCCGAGGCCAGCGGAGACCAACGCCACCTGC
 CACCACCAGGCCCGAGAGGCCAGCGGAGACCAAGGTCCCTGAGGAGATCCCCCTGAAGTGGTGCAGGA
 GTATGTGGACATCATGGAGGAGCTGCTGGGGTCTCACCTGGGGACACAGGGGAGCCTGAGGGACACCG
 GAAAAGGGCAAAGTGGAGCAGCCGAGGAAGAGGACGGGATAACCTCAGACCCGGGCCTCTGAGCTACA
 TTGACAAGCTGTGTTCCAGGAAGACTTTGTACCAAGGTGGAGGCCGTATTACCCCCGATTCTTGGA
 AGAATTGCTTTCCCAGATCCACAGATGGATTTCTTGGCCCTAAGCCAGGAGCTGGAGCAGGAGGAAGGA
 CTCACCTTGCCAGCTAGTGGAGAAGCGCTCCTGTCTTGAAGGAGAAAGGTGCGGGAGGGCAGCCC
 CTCGACATGGCACGGCCAGGTTGACTCAAGTCTTCTGAGTTTGCAGCTGGCCAAGAAGCAGCGAGAGA
 GGTCCCTGACCCCAACAAAGGGTCAGCGTGAAACCTCCCCACCCAGACAGCTGCCCAGGACCTCAG
 GGACAGGGCAGAGTGCCTACTGGCATGGCCAGGTCCGAAGACCCTGCTGTGCTTTTGGGATGTCAGGATT
 CCCCCAGGCTGAAGGCTGTCCGGCCAACCTCTCTCCCCAGGACCACAGCCACCTGCCCGGCCCTGGG
 GACCAAGGACGCCTTGGGTCTCCCTGGAGAGTCTCCTGTCAAGGAGTACATGGGCTGGCTAAGGGGTCA
 AGTGAGGAGACGGAACCTCCCTGGCATGGTCTATGTCTGGGTTCCACCACAGGCTGAGGCCCTGGAGGC
 TGTCACAGAGCCTGTCCCTTCTCGGGCTTCTCAGCCAGGAGGGAGAGGACCCAGGGAGCTTTCA
 GTCTCCATCTGCTCAGAAAAGAGGCTCAGCCCATCACCTTCTCTGCCAGCAAGTCCAAGAAGCGACCT
 CTCTTTGGAAGCCATCCCCTGCTGAAAAGACACCGCACCCAGGGCTGGGCTCAGGGTCTCTGGGGAGC
 AATCCCTGGCTTGGGGGCTGGGTGGCCCTCACAGTCTCAAAGAGAAAGGTGACCCCTTGGCCTCCAG
 GAGGAAGAAGAAGCGGCATTGTAGCCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC215918 representing NM_017561
 Red=Cloning site Green=Tags(s)


MASNGAYPVLGPGVTVPNGTSLSVFTALPFATPAPGPAHRPPLVTAVVPPAGPLVLSAFPSTPLVAGQDG
 RGP SGAGASN VFVQMRTEVGPVKPPQAQTLILTQAPLVWQAPGTL CGGVMCPPLLLLAAAPGVPVTSQV
 VGGTQACEGGWSHGLPLPPPPAAQVAPIVSPGNARPWPQGAHGEGLAPSQAKARPDDSCPKSVYENF
 RLWQHYKPLARRHLPQSPDTEALSCFLIPVLRSLARRKPTMTLEEGLWQAMREWQHTSNFDRMIFYEMAE
 KFLEFEAEEMQIQKSQWMKGPQSLPPPAPPRLEPRGPPAPEVVKQPVYLP SKDGP KAPTACLPPPRPQR
 PAETKAHLPPPRPQRPAETNAHLPPPRPQRPAETKVPEEIPPEVVQEYVDIMEELLGSHPGDTGEPEGQR
 EK GKVEQPQEEDGITS DPGLLSYIDKLCSQEDFVTKVEAVIHPRFLEELLSPDPQMDFLALSQELEQEEG
 LTLAQLVEKRLLSLKEKGCGRAAPRHGTARLDSSPSEFAAGQEAAREVPDPQQRVSVETSPQTAAQDPQ
 GQGRVRTGMARSEDPAVLLGCQDSPRLKAVRPTSPPQDHRPTCPGLGTKDALGLPGESPVKESHGLAKGS
 SEETELPGMVVYVGSHHRLRPWRLSQSPVPSSGLLSPGGRGPQALQSPSAQKRGLSPSPSPASKSKKRP
 LFGSPSPAECTHPGPGLRVSGEQLAWGLGGPSQSQRKGDPLASRRKKKRHCSQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8101_f06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_017561


ORF Size: 2268 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_017561.2</u>
RefSeq Size:	2561 bp
RefSeq ORF:	2271 bp
Locus ID:	54754
UniProt ID:	<u>A1L443</u>
Cytogenetics:	9q22.32
MW:	80.8 kDa

Product images:



Circular map for RC215918



Western blot validation of overexpression lysate (Cat# [LY413704]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC215918 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).