

## Product datasheet for **RC215922L1V**

### Spectrin beta 5 (SPTBN5) (NM\_016642) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Spectrin beta 5 (SPTBN5) (NM_016642) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Spectrin beta 5
Synonyms:	BSPECV; HUBSPECV; HUSPECV
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016642
ORF Size:	11022 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215922).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016642.3</a>
RefSeq Size:	11748 bp
RefSeq ORF:	11025 bp
Locus ID:	51332
UniProt ID:	<a href="#">Q9NRC6</a>
Cytogenetics:	15q15.1
MW:	417.2 kDa



[View online »](#)