

## Product datasheet for **RC215923L1V**

### **PACAP receptor (ADCYAP1R1) (NM\_001118) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PACAP receptor (ADCYAP1R1) (NM_001118) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PACAP receptor
Synonyms:	PAC1; PAC1R; PACAPR; PACAPRI
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001118
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215923).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001118.3</a>
RefSeq Size:	1643 bp
RefSeq ORF:	1407 bp
Locus ID:	117
UniProt ID:	<a href="#">P41586</a>
Cytogenetics:	7p14.3
Domains:	7tm_2, HormR
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

**Protein Pathways:** Neuroactive ligand-receptor interaction

**MW:** 53.31 kDa

**Gene Summary:** This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Dec 2010]