


Product datasheet for RC215929L4

CD95 (FAS) (NM_152871) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD95 (FAS) (NM_152871) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	CD95
Synonyms:	ALPS1A; APO-1; APT1; CD95; FAS1; FASTM; TNFRSF6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215929).
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	
	Cloning sites used for ORF Shuttling:
	Sgf1 ORF Mlu I GCG ATC GCC ATG NNN ACG CGT



ACCN: ORF Size: NM_152871 942 bp

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


	CD95 (FAS) (NM_152871) Human Tagged Lenti ORF Clone – RC215929L4
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Me	 2. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 152871.1</u>
RefSeq Size:	2692 bp
RefSeq ORF:	945 bp
Locus ID:	355
UniProt ID:	<u>P25445</u>
Cytogenetics:	10q23.31
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	Allograft rejection, Alzheimer's disease, Apoptosis, Autoimmune thyroid disease, Cytokine- cytokine receptor interaction, Graft-versus-host disease, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, p53 signaling pathway, Pathways in cancer, Type I diabetes mellitus
MW:	33.6 kDa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Scherkeiter CD95 (FAS) (NM_152871) Human Tagged Lenti ORF Clone – RC215929L4

Gene Summary:The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor
contains a death domain. It has been shown to play a central role in the physiological
regulation of programmed cell death, and has been implicated in the pathogenesis of various
malignancies and diseases of the immune system. The interaction of this receptor with its
ligand allows the formation of a death-inducing signaling complex that includes Fas-
associated death domain protein (FADD), caspase 8, and caspase 10. The autoproteolytic
processing of the caspases in the complex triggers a downstream caspase cascade, and leads
to apoptosis. This receptor has been also shown to activate NF-kappaB, MAPK3/ERK1, and
MAPK8/JNK, and is found to be involved in transducing the proliferating signals in normal
diploid fibroblast and T cells. Several alternatively spliced transcript variants have been
described, some of which are candidates for nonsense-mediated mRNA decay (NMD). The
isoforms lacking the transmembrane domain may negatively regulate the apoptosis mediated
by the full length isoform. [provided by RefSeq, Mar 2011]

Product images:



Circular map for RC215929L4

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US