

## Product datasheet for **RC215948L3V**

### CREM (NM\_001881) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CREM (NM_001881) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CREM
Synonyms:	CREM-2; hCREM-2; ICER
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001881
ORF Size:	411 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001881.3</a> , <a href="#">NP_001872.3</a>
RefSeq Size:	1304 bp
RefSeq ORF:	414 bp
Locus ID:	1390
UniProt ID:	<a href="#">Q03060</a>
Cytogenetics:	10p11.21
Domains:	pKID, BRLZ
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 14.9 kDa

**Gene Summary:** This gene encodes a bZIP transcription factor that binds to the cAMP responsive element found in many viral and cellular promoters. It is an important component of cAMP-mediated signal transduction during the spermatogenic cycle, as well as other complex processes. Alternative promoter and translation initiation site usage allows this gene to exert spatial and temporal specificity to cAMP responsiveness. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene, with some of them functioning as activators and some as repressors of transcription. [provided by RefSeq, Jul 2008]