

Product datasheet for **RC215951A1V**

Human FBXW9 (NM_032301) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FBXW9 (NM_032301) AAV Particle
Tag:	Myc-DDK
Symbol:	FBXW9
Synonyms:	Fbw9; MEC-15
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC215951 representing NM_032301
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGCTTCCCCTAGGGCGGTGCGATGATTCCCGCACCTGGGACGATGACTCGGACCCAGAGTCAGAGA
 CAGACCCAGACGCGCAGGCCAAGGCCTACGTGGCCCGCGTTCTCAGTCCGCCAAAATCCGGGCTGGCGTT
 CTCGCGCCCTCGCAGCTATCCACACCCGCCGCGTCCCGAGCGCTTCGGAGCCTCGGGCCGCGTCCAGG
 GTTTCGGCCGTAAGTGAGCCGGGCTTCTGAGCCTTCCCCGGAGCTGCTGCTCGAGATCTGCTCCTACC
 TGGACGCCCGCTCGTGTCCACGTCTGTGCGGGGTGTGCCACGCGCTCCGCGACCTCGTGTCTGACCA
 TGTACCTGGAGGCTACGCGCGTACGCCGCGTACGCGCGCCCTACCCAGTGGTGAAGAGAAGAATTT
 GACTGGCCGGCAGCCTGCATTGCGCTGGAGCAGCACCTGTCCGCTGGGAGAGGATGGGCGTGGTCTG
 AATACTTCTGCTGGCCGAAGGCCACGTGGCTTCCGTTGACTCAGTGTCTGCTCCAGGTGGTCACT
 CTGTCTGTGCGGCTCCCGAGATCGAACGTCAACTTGTGGGACCTGCGGACGCTGGGACGGAGTCCAAC
 CAGGTTCTGATCAAGACCTTAGGCACTAAGCGAAATAGTACCCATGAGGGCTGGGTGTGGTCACTGGCAG
 CGCAGGACCACCGCTGTGCTCCGGCTCCTGGGACAGCACAGTGAAGCTCTGGGACATGGCAGCGGATGG
 GCAGCAGTTCGGCGAGATAAAGGCCAGCTCAGCCGTGCTGTGCCTCTCTACCTGCCTGACATCCTGGTG
 ACTGGCACCTATGACAAGAAGGTGACCATCTACGACCCAGAGCCGGCCAGCCCTGTTGAAGCACCAGC
 AACTACACTCCAGACCCGTGCTGACCTGTGGCGGATGACCGGCACATCATCTCAGGCAGCGAGGACCA
 CACCCTGGTGGTGGTGGACGCCGAGCCAACAGCGTCTGCAGCGTCTGCAGTGGACTCTACCTGCTC
 TGCATGCTCTACCAGGAACCCAGCTCTGGGCTGGTGACAACAGGGCCTGCTGCAGCTCTTCGCCAACC
 GCAACGGCTGCTTCCAGCTTATCCGGTCTTTGATGTGGGCCACAGCTTCCCATCACTGGGATCCAGTA
 CTCGCTGGGAGCCTTGTACACCACATCCACTGACAAGACCATCCGGGTGCAGTGGCCACAGACCCACCA
 AGGACCATTTGACCCGAAGGCATGACATGGGCTCAATAGGGTCTGTGCTGAGGGCAACCTGGTGGTGG
 CCGGCTCTGGAGACCTGTGCTAGAGGTCTGGAGGCTGCAGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC215951 representing NM_032301
 Red=Cloning site Green=Tags(s)

MELPLGRCDSDRTWDDSDPESETDPDAQAKAYVARVLSPPKSGLAFSRPSQLSTPAASPSASEPRAASR
 VSAVSEPLLSLPPELLLEICSYLDARLVHLVSRVCHALRDLVSDHVTWRLRALRRVRAPYPVVEEKNF
 DWPAACIALEQHLRWAEDGRWEYFCLAEGHVASVDSVLLQGGSLCLSGSRDRNVNLWDLRQLGTESN
 QVLIKTLGTRNSTHEGWVWSLAAQDHRVCSGWDSTVKLWDMAADGQQFGEIKASSAVLCLSYLPDILV
 TGTYDKKVITYDPRAGPALLKHQQLHSRPVLTLLADDRHIIISGSEDHTLVVVDRRANSVLQRLQLDSYLL
 CMSYQEPQLWAGDNQGLLHVAFANRNGCFQLIRSFVGHSPITIGIQSVGALYTTSTDKTIRVHVPTDPP
 RTICTRRHDNGLNRVCAEGNLVVAGSGDLSLEVWRLQA

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species:

Human

Serotype:

AAV-2

ACCN:

NM_032301

ORF Size:

1374 bp

Buffer:

PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032301.2</u>
RefSeq Size:	1750 bp
RefSeq ORF:	1377 bp
Locus ID:	84261
UniProt ID:	<u>Q5XUX1</u>
Cytogenetics:	19p13.13
MW:	50.6 kDa