

## Product datasheet for **RC215988L4V**

### **FAM111B (NM\_198947) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FAM111B (NM_198947) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM111B
Synonyms:	CANP; POIKTMP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198947
ORF Size:	2202 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215988).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198947.1</a>
RefSeq Size:	3574 bp
RefSeq ORF:	2205 bp
Locus ID:	374393
UniProt ID:	<a href="#">Q6SJ93</a>
Cytogenetics:	11q12.1
MW:	84.7 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein with a trypsin-like cysteine/serine peptidase domain in the C-terminus. Mutations in this gene are associated with an autosomal dominant form of hereditary fibrosing poikiloderma (HFP). Affected individuals display mottled pigmentation, telangiectasia, epidermal atrophy, tendon contractures, and progressive pulmonary fibrosis. Alternative splicing results in multiple transcript variants encoding distinct isoforms. A paralog of this gene which also has a trypsin-like peptidase domain, FAM111A, is located only 16 kb from this gene on human chromosome 11q12.1. [provided by RefSeq, Apr 2014]