

## Product datasheet for **RC215991L3V**

### CCDC91 (NM\_018318) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC91 (NM_018318) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC91
Synonyms:	HSD8; p56
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018318
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215991).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018318.3</a>
RefSeq Size:	2352 bp
RefSeq ORF:	1326 bp
Locus ID:	55297
UniProt ID:	<a href="#">Q7Z6B0</a>
Cytogenetics:	12p11.22
MW:	49.8 kDa



[View online »](#)

**Gene Summary:**

Involved in the regulation of membrane traffic through the trans-Golgi network (TGN). Functions in close cooperation with the GGAs in the sorting of hydrolases to lysosomes. [UniProtKB/Swiss-Prot Function]