

## Product datasheet for **RC215998L4V**

### **SYTL5 (NM\_138780) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SYTL5 (NM_138780) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SYTL5
Synonyms:	slp5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138780
ORF Size:	2190 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215998).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138780.1</a>
RefSeq Size:	2193 bp
RefSeq ORF:	2193 bp
Locus ID:	94122
UniProt ID:	<a href="#">Q8TDW5</a>
Cytogenetics:	Xp11.4
MW:	81.3 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the synaptotagmin-like (Slp) protein family, which contains a unique homology domain at the N-terminus, referred to as the Slp homology domain (SHD). The SHD functions as a binding site for Rab27A, which plays a role in protein transport. Expression of this gene is restricted to placenta and liver, suggesting that it might be involved in Rab27A-dependent membrane trafficking in specific tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]